





Thermarmour Emergency Poncho

- Material:** Non-conductive bi-laminate foil & needle punched polyester fleece
- Colour:** Blue outer / White inner
- Size:** 100cm x 137cm (can be cut to size)
- Product code:** OB210LGE-----TAEM55
- Washable:** No – single use
- Purchase online:** www.interweavetextiles.com/products/patient-warming



Thermarmour is a trading brand of Interweave Textiles Ltd.
Old Power Way, Lowfields Business Park, Elland, West Yorkshire, UK, HX5 9DE, Tel: 01422 372333

Product features


No heat source required
One size fits most
Shower proof
High visibility levels
Latex-free
Recyclable
Standalone single use
Can be cut to size

Ideal for


Emergency services
Search & rescue teams
First responders
Military personnel
Private medical providers
Recovery vehicles and breakdowns
Outdoor enthusiasts



1 Using the provided slit, tear open the packaging, pull out poncho and unfold.



2 Place over the patient's head, foil on the outside.



3 Can be cut to fit the patient [* see note]

Thermarmour is a trading brand of Interweave Textiles Ltd.
Old Power Way, Lowfields Business Park, Elland, West Yorkshire, UK, HX5 9DE Tel: 01422 372333