



Thermarmour Emergency Baby Wrap

Clinically proven and conforms to World Health Organization (WHO) “chain warming” guidelines.

- Material:** Non-conductive bi-laminate foil & needle punched polyester fleece
- Colour:** Blue outer / White inner
- Size:** 55cm x 57cm [can be cut to size]
- Product code:** 0B202-55-57---TAEM55
- Washable:** No – single use
- Buy online:** www.interweavetextiles.com/products/patient-warming/



Thermarmour is a trading brand of Interweave Textiles Ltd.
Old Power Way, Lowfields Business Park, Elland, West Yorkshire, UK, HX5 9DE, Tel: 01422 372333

Product features

Suitable for full-term babies
Helps maintain thermoregulation
Can be pre-warmed in blanket warmer
Waterproof & wind proof
X-ray compatible
Latex-free
Recyclable
Standalone single use
Can be cut to size

Ideal for

Emergency services
Search & rescue teams
First responders
Private medical providers
Paediatrics wards


Instructions:



1) After opening packaging, place baby in the middle of the Baby Wrap



2) Wrap baby in the Baby Wrap and secure



Thermarmour is a trading brand of Interweave Textiles Ltd.
Old Power Way, Lowfields Business Park, Elland, West Yorkshire, UK, HX5 9DE Tel: 01422 372333