



## Thermarmour Emergency Blanket 2m x 1.5m

**Material:** Non-conductive bi-laminate foil & needle punched polyester fleece


**Colour:** Blue outer / White inner

**Size:** 2m x 1.5m

**Product code:** OB000150200---TAEM55

**Washable:** No – single use

**Purchase online:** [www.interweavetextiles.com/products/patient-warming](http://www.interweavetextiles.com/products/patient-warming)



Interweave Textiles Ltd, Old Power Way, Lowfields Business Park, Elland, West Yorkshire, UK, HX5 9DE, Tel: 01422 372333

Document valid as at 04/09/2023


Page 1 of 2 - Revision no. 1

## Product features


Helps maintain thermoregulation  
Shower proof  
High visibility levels  
No heat source required  
Latex-free  
Recyclable  
Standalone single use  
Can be cut to size

## Ideal for

Emergency services  
Search & rescue teams  
First responders  
Military personnel  
Private medical providers  
Recovery vehicles  
Breakdown assistance



- 1.** Using the provided slit, tear open the packaging and pull out the blanket. Unfold to full size.



- 2.** Wrap the patient in the blanket with the foil side on the outside