



## Thermarmour Medical Blanket 1m x 1.5m

**Material:** Non-conductive bi-laminate foil & needle punched polyester fleece


**Colour:** Light blue outer / White inner

**Size:** 1m x 1.5m

**Product code:** OB000100150---TAMM55

**Washable:** No – single use

**Purchase online:** [www.interweavetextiles.com/products/patient-warming](http://www.interweavetextiles.com/products/patient-warming)



Thermarmour is a trading brand of Interweave Textiles Ltd.  
Old Power Way, Lowfields Business Park, Elland, West Yorkshire, UK, HX5 9DE, Tel: 01422 372333


## Product features

No heat source required  
X-ray compatible  
Latex-free  
Recyclable  
Standalone single use  
Can be cut to size

## Ideal for

Surgery: Peri, intra and post operatively  
Suitable for use with diathermy & cautery treatments  
Situations where forced-air warming is not applicable or cost-effective

## Instructions



- 1] Using the provided slit, tear open the packaging, pull out the blanket & unfold.
- 2] Either cover or wrap the patient completely with the foil on the outside, 30 minutes prior to any procedure.
- 3] Can be cut to fit the patient as long as the majority of the thermal pockets in direct contact with the patient remain intact.

Thermarmour is a trading brand of Interweave Textiles Ltd.  
Old Power Way, Lowfields Business Park, Elland, West Yorkshire, UK, HX5 9DE Tel: 01422 372333