

# THE AGILE™ HYPOTHERMIA EVACUATION BAG

COMES WITH MARCH CHECKLIST™, PROTECTIVE LINER AND MULTIFUNCTIONAL GEAR BAG

## ✓ THE AGILE SYSTEM

- Modular system, informed by military operational experience and the reality of casualty care
- Advanced passive warming to aid hypothermia prevention
- Lightweight, robust, compact - pack size matches a medical backpack
- Reusable and easy to clean
- Engineered using reliable, field-proven components from leaders in outdoor and military equipment
- THERMARMOUR® range continuously improving and adapting for military and emergency use
- Manufactured in the UK

## ✓ Easy handling under pressure in harsh environments (MASCAL)

- Evacuation bag - rolled up, ready for use, compressed with Velcro strap
- Can be shaken out of the gear bag with one hand and is ready for use in seconds
- Large icon at the head end for orientation
- MARCH Checklist™ on the chest area weatherproof, writable, removable for documentation and handover [CRM]
- Gear bag secured against loss with paracord, suitable as a dump pouch / CAS equipment bag

## ✓ THE EVACUATION BAG

- **Material:** Patented multilaminate THERMARMOUR® Xtreme to minimise the effect of hypothermia
- **Colour:** Outside olive-green; inside white, making it easier to detect blood /fluid loss
- **Total size unrolled:** 235 cm x 106 cm **Total size rolled:** 35 cm x 20 cm **Weight:** 2.5 kg  
Suitable for use with Therm-a-Rest™ or rollmat
- **Water & wind-resistant shell** for effective environmental protection and **wipeable inner surface**
- **Casualty access:** Front can be entirely opened for quick/easy loading/unloading CAS and medical treatment. All around heat and moisture-tight zipper for easy loading and treatment and secure transport of CAS. Features hood with elastic hem to seal tightly around the face
- **Colour coding:** 4 pairs of colour-coded zippers allow treatment openings wherever needed [SEAL / NDC / DRAIN / I.V. / I.O. / ACTIVE WARMING]
- **Patient handling:** 8 carrying handles for safe and efficient evacuation - also suitable for 2 rescuers. Tubular Webbing [25 mm, 20 kN, EN 565] cross hatched across the entire floor area for stabilised transport of the CAS [CQB, handling on Platform CASEVAC MEDEVAC]
- **Compression strap:** Velcro compression strap for ease of packing and transportation
- **Classification:** MDR [EU]2017/745 Medical device class I

# THE AGILE

## ✓ THE AGILE MARCH CHECKLIST™

- **MARCH FIRST**
- Specifically designed MARCH Checklist™ to ensure structured treatment routine, documentation and transfer: Clear, weatherproof, detachable
- Fixed with velcro at chest height
- Hole rivet for pen attachment

## ✓ PROTECTIVE LINER

- Replaceable
- Highly absorbent
- Dimensions 70 cm x 100 cm, fitted to base
- Fast to change for the next CAS

## ✓ MULTIFUNCTIONAL GEAR BAG

- Robust Cordura according to TL: Tactical IR equipment. Water, oil, and dirt repellent
- **Length:** Adaptable 40 - 75 cm **Diameter:** 25 cm
- Volume adjustable: 8 litres (evacuation bag only) or 14 litres (for additional gear)
- Pack size when compressed with evacuation bag, including protective liner 25 cm x 40 cm
- Roll top: waterproof closure, quick opening and closing
- Adjustable, detachable shoulder strap
- Can be attached to the med bag or used as a rapid access trauma pouch

## ✓ ADDITIONAL EQUIPMENT

### Multifunctional Poncho

- Additional THERMARMOUR® hypothermia prevention. Can be used for weather protection, as an isolation pad for treatment on the ground, and as a thermal cover for limbs during medical treatment.

### Replacement sets

- MARCH Checklist™, 5 pcs/set
- Protective liner, 5 pcs/set

SCAN THE QR CODE TO VIEW PRODUCTS

☎ 01422 411463

✉ sales@thermarmour.com



www.thermarmour.com